




SMART CRYO-TOOLS

Zeus Smart Cryo Tools takes your freezing process to the next level.



SMART CRYO-TOOLS

Zeus Smart Cryo Tools takes your freezing process to the next level.

The octagonal goblet design allows for the radial arrangement of **8 straws**, saving space.




Thanks to its ergonomic straw body and soft, transparent, and ideally structured straw tip, you can load oocyte or embryo microdroplets onto the straw under a stereo microscope and perform a safe freezing process.



The radial arrangement provides a significant advantage and added reliability by enabling visual inspection of the straw labels at the tank before dissolving and, more importantly, confirming the barcode with a wireless reader.

Smart straws, with their adjustable lengths and triangular design can **increase the existing volume of your tank by up to 4 times.**




ZEUS
Cryo-Tools Cane



Traditional
Cane


Zeus Smart Tools canes fit alongside existing canes in the same canister, instantly increasing tank capacity without the need for a new liquid nitrogen tank.

Confirm embryos at every freezing or thawing step using a wireless barcode reader.



During thawing, shortened ZEUS smart straws are easily removed with an extender grabber pen.

ZEUS Cryotools maximizes safety with real-time barcode verification of each straw.



For PGT-tested IVF cycles, correct embryo identification during thawing is as critical as patient confirmation.