



**B-TEST**  
Testing For Life

**BIÇAKCILAR**  
**LAB**

Osmangazi Mah. Gazi Cad. No:21 Esenyurt İstanbul  
Tel: +90 (212) 689 02 20 Fax: +90 (212) 689 02 29  
labinfo@bicakcilar.com

|                                       |            |                                                                |                                    |
|---------------------------------------|------------|----------------------------------------------------------------|------------------------------------|
| <b>Request-Tender no</b>              | TLF-21-947 | <b>STERILITY TEST REPORT</b><br><b>(European Pharmacopeia)</b> | <br>AB-0052-T<br>M21-1716<br>07-21 |
| <b>Receiving Date</b>                 | 22.06.2021 |                                                                |                                    |
| <b>Test Initiation Date</b>           | 24.06.2021 |                                                                |                                    |
| <b>Test Final Date</b>                | 08.07.2021 |                                                                |                                    |
| <b>Publication Date of the Report</b> | 08.07.2021 |                                                                |                                    |

#### **CUSTOMER-ADDRESS**

ZEUS BIOSCIENCE BİYOMEDİKAL A.Ş.  
ÇANKAYA MAH. CİNNAH CADDESİ 47/2 06690 ÇANKAYA/ANKARA

#### **SAMPLE TRANSMISSION TYPE**

It was sent by the customer.

#### **METHOD DESCRIPTION**

Sterility test was performed under Laminar Air Flow by direct transfer method.

Sample or sample prepared as per SIP was transferred to the jar, filled with TSB growth medium enough to cover the sample.

Prepared, TSB medium with product was incubated at (20-25)°C.

Other were transferred to the jar, filled with TFM growth medium enough to cover the whole product.

Prepared, TFM medium was incubated at (30-35)°C.

The jars containing the products were controlled for growth of microorganisms for 14 days.

#### **REFERENCE STANDARD**

European Pharmacopoeia (EP 2.6.1.)

#### **SAMPLE IDENTIFICATION - SAMPLE NO**

ZEUS Electroliquid Solution /LotNo:100006 /RefNo:ZPS19-003 / M21-1716-1



**B-TEST**  
Testing For Life

**BIÇAKCILAR**  
**LAB**

Osmangazi Mah. Gazi Cad. No:21 Esenyurt İstanbul  
Tel: +90 (212) 689 02 20 Fax: +90 (212) 689 02 29  
labinfo@bicakcilar.com

AB-0052-T

M21-1716

07-21

#### ▮ SAMPLE QUANTITY

2 Piece

#### ▮ SAMPLE ITEM PORTION (SIP)

NA

#### ▮ TEST RESULTS

|                                            | for TSB | for TFM |
|--------------------------------------------|---------|---------|
| Number of product controlled               | 1       | 1       |
| Number and date of product growth observed | 0       | 0       |

It has been declared by the customer that the product does not have reproductive preventive or supportive

#### ▮ DECLARATION OF CONFORMITY

Declaration of conformity is not provided according to the test method.

#### ▮ LIMITS

Shall be determined by the customer

#### ▮ DEVIATIONS

NA

#### ▮ REVISION HISTORY

| Rev No - Date   | Revision Reason |
|-----------------|-----------------|
| 00 - 08.07.2021 | New Publication |



**B-TEST**  
Testing For Life

**BIÇAKCILAR**  
**LAB**

Osmangazi Mah. Gazi Cad. No:21 Esenyurt İstanbul  
Tel: +90 (212) 689 02 20 Fax: +90 (212) 689 02 29  
labinfo@bicakcilar.com

AB-0052-T

M21-1716

07-21

## II APPROVALS

|                      | Name - Title                  | Signature - Date                         |
|----------------------|-------------------------------|------------------------------------------|
| Laboratory Approval: | Ayşegül İSFENDİYAR<br>CHEMIST |                                          |
| Quality Approval:    | Özenç EFE ÖZTÜRK<br>CHEMIST   |                                          |
| Seal                 |                               | <b>BIÇAKCILAR</b><br><b>LABORATUVARI</b> |

Test results are valid only for the tested samples identified in this test report

The laboratory is not responsible for deviations in analysis results if the declared information is incomplete or inaccurate, or if the customer is requested to work it although the samples are not delivered in accordance with the acceptance criteria.

This report shall not be reproduced other than in full except with the permission of Bıçakcılar Laboratory.

Testing reports without signature and seal are not valid.

This document is signed with a secure electronic signature according to the electronic signature law numbered 5070.